



Medical mounting arms and mounting systems

PREMIUM QUALITY WITH
EXCEPTIONAL FLEXIBILITY



TECHNICAL DETAILS

rm-port	single	double	double
Lengths	250 mm / 325 mm	250 mm + 250 mm	325 mm + 325 mm
Monitor weight	up to 18 kg / 39.6 lbs	up to 18 kg / 39.6 lbs	up to 18 kg / 39.6 lbs
Load capacity	up to 23 kg / 50.6 lbs	up to 23 kg / 50.6 lbs	up to 18 kg / 39.6 lbs

mf-port	single	double
Lengths	350 mm	350 mm + 350 mm
Monitor weight	up to 18 kg / 39.6 lbs	up to 18 kg / 39.6 lbs
Load capacity	up to 23 kg / 50.6 lbs	up to 23 kg / 50.6 lbs

lf-port	support arm	support arms on pivot arm
Lengths	225 mm	225 mm + 250 mm 225 mm + 325 mm
Monitor weight	up to 7 kg / 15.4 lbs	up to 8 kg / 17.6 lbs up to 8 kg / 17.6 lbs
Load capacity	up to 7 kg / 15.4 lbs	up to 18 kg / 39.6 lbs up to 18 kg / 39.6 lbs

flexion-port	support arms	support arms on pivot arm
Lengths	300 mm / 400 mm	300 mm + 350 mm 400 mm + 350 mm
Monitor weight	up to 20 kg / 44 lbs	up to 14 kg / 30.8 lbs up to 14 kg / 30.8 lbs
Load capacity	up to 20 kg / 44 lbs	up to 14 kg / 30.8 lbs up to 14 kg / 30.8 lbs

Get in touch »

Europe: +49 89 614425-0 | sales@itd-cart.com
America: +1 800 947-3901 | salesusa@itd-cart.com
China: +86 21 5031 8818 | saleschina@itd-cart.com
India: +91 86 0015 4064 | salesindia@itd-cart.com
Australia: +61 3 9330 5050 | salesaustralia@itd-cart.com

www.itd-cart.com




MOUNTING SYSTEMS

Medical mounting system modul-port



Different iTD profiles –
in lengths from 320 mm to 1920 mm



A wide range of monitor mounts –
to accommodate standard monitors



Large portfolio of accessories and components –
for various applications



DS.0292.799 – 10/2026. Printed in Germany. All rights to make changes reserved.

OVERVIEW MEDICAL MOUNTING ARMS

Pivot arms rm-port

Can be swiveled, tilted and rotated –
wide range of variants with horizontal
swivel range of $\pm 90^\circ$, tilt range of up to
up to $\pm 20^\circ$ and rotatable by $\pm 5^\circ$


Various arm lengths –
swivel arms single in 250 mm and 325 mm
and swivel arms double in 250 mm + 250 mm
and 325 mm + 325 mm

Effortless swiveling of even heavy screens –
high ease of use and reliability
due to maintenance-free joints

Numerous interfaces –
for mounting standard monitors

Numerous interfaces –
for mounting on the hospital
infrastructure

Integrated mains lead cleats –
protect the cable against being pulled out



Pivot arms mf-port

Various lengths –
available as single or double arm,
customizable in length (L1 and L2)
and number of joints

Adjustable swivel joints –
adjustment of friction within the joints


Numerous interfaces –
for mounting standard monitors
like VESA 75 / 100, universal adapter,
table top mount adapter

Numerous interfaces –
for mounting on the hospital
infrastructure

CUSTOMIZED SOLUTIONS



Variable height support arms lf-port



Different lengths –
enable a flexible and ergonomic
design of the workplace

High level of operating comfort and reliability –
due to integrated spring force support

Large-volume cable ducts –
for data and power cables

Numerous interfaces –
for mounting on the hospital
infrastructure

Rapid assembly –
specially developed rapid fasteners

Variable height support arms flexion-port

Different lengths –
enable a flexible and ergonomic
design of the workplace

Effortless lifting of heavy monitors –
high operating comfort and reliability
due to spring force assistance

Numerous interfaces –
for mounting standard monitors

Large-volume cable ducts –
for data and power cables

Numerous interfaces –
for mounting on the hospital
infrastructure

