

classic-cart Profi

Planning Support

Equipment cart

Power columns: yes GW. _____

On and off switch yes no

Castor conversion kit: ZV. _____

Castor diameter: D _____ mm

Conductivity: yes no

Isolating transformers: 660 VA with 4-way
EURO socket strip
 1000 VA with 8-way
EURO socket strip
 1600 VA with 8-way
EURO socket strip

Insulation monitor yes no

Socket strip: ZV. _____

Rear panel: yes no

Accessories:


Design

Base colour: RAL 7035 (Light grey)

Trim colour: .901 (RAL 5021 Water blue)
 .902 (RAL 5013 Cobalt blue)
 .903 (RAL 7035 Light grey)

Please note:

Just draw in your selection of accessories at your choice of position between the two dotted lines. Accessories like shelves, drawers etc. for **classic-cart** can be assembled starting **at the second unit (at 64 mm)** from the base.



Dimensions

Total height TG:

Units (U), 1 U = 32 mm	Castors D 100 mm	Castors D 125 mm
21 U	829 mm	857 mm
30 U	1117 mm	1145 mm
40 U	1437 mm	1465 mm

Installation width: 500 / 600 mm

Total width: 614 / 714 mm

Total depth: 625 mm

Project _____

Company _____

Street _____

City _____

Fon _____

Fax _____

E-Mail _____

Processor _____

iTD GmbH
Lilienthalstr. 19
85579 Neubiberg
Tel. +49 89 61 44 25 100
E-Mail: sales@itd-cart.com
Web: www.itd-cart.com

