

pro-cart

Planning Support

Equipment

Castor diameter: D _____ mm

Conductivity: yes no

Isolating transformers: 660 VA with 4-way EURO socket strip
 1200 VA with 6-way EURO socket strip

Insulating monitor: yes no


Socket strip: ZV. _____

Accessories: _____

Design

Base colour: RAL 7035 (Light grey)

Trim colour: .903 (RAL 7035 Light grey)
 .908 (RAL 5015 Sky blue)
 .909 (RAL 1003 Signal yellow)



Dimensions

Total height TH:

Units (U), 1 U = 32 mm	Castors D 100 mm	Castors D 125 mm
21 U	884 mm	909 mm
30 U	1172 mm	1197 mm
35 U	1332 mm	1357 mm

Total width: 610 mm

Total depth: 620 mm

Project _____

Company _____

Street _____

City _____

Fon _____

Fax _____

E-Mail _____

Processor _____

ITD GmbH
 Lilienthalstr. 19
 85579 Neubiberg, Germany
 Tel. +49 89 61 44 25 100
 Fax +49 89 61 44 25 200
 E-Mail: sales@itd-cart.com
 Web: www.itd-cart.com



Please note:

Just draw in your selection of accessories at your choice of position between the two dotted lines. Accessories like shelves, drawers etc. for **pro-cart** can be assembled starting **at the first unit (at 32 mm)** from the base.