

vexio-cart Plus

Planning Support

Equipment

Castor diameter: D _____ mm

Conductivity: yes no

Isolating transformers: 660 VA with 4-way EURO socket strip
 1200 VA with 6-way EURO socket strip

Insulating monitor: yes no


Socket strip: ZV. _____

Accessories: _____

Design

Base colour: RAL 7035 (Light grey)

Trim colour: .903 (RAL 7035 Light grey)
 .908 (RAL 5015 Sky blue)
 .909 (RAL 1003 Signal yellow)



35 U (1120 mm)

30 U (960 mm)

25 U (800 mm)

21 U (672 mm)

15 U (480 mm)

10 U (320 mm)

5 U (160 mm)

1 U (32 mm)

Dimensions

Total height TH:

| Units (U), 1 U = 32 mm | Castors D 100 mm | Castors D 125 mm |
|---------------------------|------------------|------------------|
| 21 U | 884 mm | 909 mm |
| 30 U | 1172 mm | 1197 mm |
| 35 U | 1332 mm | 1357 mm |

Total width: 610 mm

Total depth: 620 mm

Project _____

Company _____

Street _____

City _____

Fon _____

Fax _____

E-Mail _____

Processor _____

Please note:

Just draw in your selection of accessories at your choice of position between the two dotted lines. Accessories like shelves, drawers etc. for **pro-cart** can be assembled starting **at the first unit (at 32 mm)** from the base.

ITD GmbH
 Lilienthalstr. 19
 85579 Neubiberg, Germany
 Tel. +49 89 61 44 25 100
 Fax +49 89 61 44 25 200
 E-Mail: sales@itd-cart.com
 Web: www.itd-cart.com

