

# classic-cart Profi

## Planning Support

### Equipment

Power columns:  yes GW. \_\_\_\_\_

On and off switch  yes  no

Castor conversion kit: ZV. \_\_\_\_\_

Castor diameter: D \_\_\_\_\_ mm

Conductivity:  yes  no

Isolating transformers:  660 VA with 4-way EURO socket strip  
 1000 VA with 8-way EURO socket strip  
 1600 VA with 8-way EURO socket strip

Insulating monitor  yes  no

Socket strip: ZV. \_\_\_\_\_

Rear panel:  yes  no

Accessories: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_


### Design

Base colour: RAL 7035 (Light grey)

Trim colour:  .901 (RAL 5021 Water blue)  
 .902 (RAL 5013 Cobalt blue)  
 .903 (RAL 7035 Light grey)  
 \_\_\_\_\_

### Please note:

Just draw in your selection of accessories at your choice of position between the two dotted lines. Accessories like shelves, drawers etc. for **classic-cart** can be assembled starting **at the second unit (at 64 mm)** from the base.



## Dimensions

### Total height TG:

Units (U), 1 U = 32 mm	Castors D 100 mm	Castors D 125 mm
21 U	829 mm	857 mm
30 U	1117 mm	1145 mm
40 U	1437 mm	1465 mm

Installation width: 500 / 600 mm

Total width: 614 / 714 mm

Total depth: 625 mm

Project \_\_\_\_\_

Company \_\_\_\_\_

Street \_\_\_\_\_

City \_\_\_\_\_

Fon \_\_\_\_\_

Fax \_\_\_\_\_

E-Mail \_\_\_\_\_

Processor \_\_\_\_\_

ITD GmbH  
 Lilienthalstr. 19  
 85579 Neubiberg

Tel. +49 89 61 44 25 100  
 Fax +49 89 61 44 25 200

E-Mail: sales@itd-cart.com  
 Web: www.itd-cart.com

