

compact-cart Profi

Planning Support

Equipment

Power columns: Standard

On and off switch: yes no

Power-Unit: ja nein

Castor conversion kit: ZV. _____

Castor diameter: D _____ mm

Conductivity: yes no

Trenntransformator: 1200 VA with 6-way EURO socket strip
 2000 VA with 8-way EURO socket strip

Insulation monitor: yes no

Socket strip: ZV. _____

Rear doors: yes no

Accessories: _____


Design

Base colour: RAL 7035 (Light grey)

Trim colour: .903 (RAL 7035 Light grey)
 .908 (RAL 5015 Sky blue)
 .909 (RAL 1003 Signal yellow)

Please note:

Just draw in your selection of accessories at your choice of position between the two dotted lines. Accessories like shelves for **compact-cart** can be assembled starting **at the second unit (at 64 mm)** from the base, drawers **at the first unit (at 32 mm)**.



Dimensions

Total height TH:		
Units (U), 1 U = 32 mm	Castors D 100 mm	Castors D 125 mm
30 U	1146 mm	1171 mm
35 U	1306 mm	1331 mm
40 U	1466 mm	1491 mm

Installation width:	450 / 500 / 650 mm
Total width:	618 / 668 / 818 mm
Total depth:	624 mm

Project _____

Company _____

Street _____

City _____

Fon _____

Fax _____

E-Mail _____

Processor _____

ITD GmbH
 Lilienthalstr. 19
 85579 Neubiberg, Germany
 Tel. +49 89 61 44 25 100
 Fax +49 89 61 44 25 200
 E-Mail: sales@itd-cart.com
 Web: www.itd-cart.com

