

uni-cart


with 4 stabilizers Planning Support

Equipment

Castor diameter:	D	mm
Conductivity:	<input type="checkbox"/> yes	<input type="checkbox"/> no
Isolating transformer:	<input type="checkbox"/> 660 VA with EURO connection cable	
Socket strip:	ZV.	_____
Counterweight kit:	<input type="checkbox"/> RS.0344.999	
Accessories:	_____ _____ _____	

Design

Basic colour:	RAL 7035 (Light grey)	
Trim colour:	<input type="checkbox"/> .901 (RAL 5021 Water blue)	
	<input type="checkbox"/> .902 (RAL 5013 Cobalt blue)	
	<input type="checkbox"/> .903 (RAL 7035 Light grey)	
	<input type="checkbox"/> _____	



Dimensions

Total height TH:

Units (U), 1 U = 32 mm	Castors D 75 mm	Castors D 100 mm
21 U	820 mm	857 mm
30 U	---	1145 mm
35 U	---	1305 mm

Total width: 540 mm

Total depth: 580 mm

Project _____

Company _____

Street _____

City _____

Fon _____

Fax _____

E-Mail _____

Processor _____

Please note:

Just draw in your selection of accessories at your choice of position between the two dotted lines. A bracket is always mounted underneath the shelves or drawers in **uni-cart**. This means that the storage platform is only usable with the **third unit (96 mm)** or in the case of **drawer with the sixth unit (192 mm)** from the bottom.

ITD GmbH
Lilienthalstr. 19
85579 Neubiberg GERMANY
Tel. +49 89 61 44 25 100 Fax
+49 89 61 44 25 200
E-Mail: sales@itd-cart.com
Web: www.itd-cart.com

