

uni-cart

with 5 stabilizers

Planning Support

Equipment

Castor diameter: D _____ mm

Conductivity: yes no

Socket strip: ZV. _____

Counterweight kit: RS.0349.999


Accessories: _____

Please note: No shelves or drawers may be mounted on the **uni-cart** with 5 stabilizers.

Design

Basic colour: RAL 7035 (Light grey)

Trim colour: .901 (RAL 5021 Water blue)
 .902 (RAL 5013 Cobalt blue)
 .903 (RAL 7035 Light grey)



Dimensions

Total height TH:

| Units (U), 1 U = 32 mm | Castors D 75 mm | Castors D 100 mm |
|---------------------------|-----------------|------------------|
| 21 U | 820 mm | 857 mm |
| 30 U | --- | 1145 mm |
| 35 U | --- | 1305 mm |

Total width: 670 mm

Total depth: 650 mm

Project _____

Company _____

Street _____

City _____

Fon _____

Fax _____

E-Mail _____

Processor _____

Please note:

Just draw in your selection of accessories at your choice of position between the two dotted lines. No shelves or drawers may be mounted on the **uni-cart** with 5 stabilizers. You should use **uni-cart** with 4 stabilizers for this configuration.

ITD GmbH
 Lilienthalstr. 19
 85579 Neubiberg, GERMANY
 Tel. +49 89 61 44 25 100
 Fax +49 89 61 44 25 200
 E-Mail: sales@itd-cart.com
 Web: www.itd-cart.com

