

# vexio-cart

## Planning Support

### Equipment

**Castor diameter:** D \_\_\_\_\_ mm

**Conductivity:**  yes  no

**Isolating transformers:**  660 VA with 4-way EURO socket strip  
 1200 VA with 6-way EURO socket strip

**Insulating monitor:**  yes  no


**Socket strip:** ZV. \_\_\_\_\_

**Accessories:** \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

### Design

**Base colour:** RAL 7035 (Light grey)

**Trim colour:**  .903 (RAL 7035 Light grey)  
 .908 (RAL 5015 Sky blue)  
 .909 (RAL 1003 Signal yellow)  
 \_\_\_\_\_



### Dimensions

**Total height TH:**

Units (U), 1 U = 32 mm	Castors D 100 mm	Castors D 125 mm
21 U	884 mm	909 mm
30 U	1172 mm	1197 mm
35 U	1332 mm	1357 mm

**Total width:** 610 mm

**Total depth:** 620 mm

Project \_\_\_\_\_

Company \_\_\_\_\_

Street \_\_\_\_\_

City \_\_\_\_\_

Fon \_\_\_\_\_

Fax \_\_\_\_\_

E-Mail \_\_\_\_\_

Processor \_\_\_\_\_

ITD GmbH  
 Lilienthalstr. 19  
 85579 Neubiberg, Germany  
 Tel. +49 89 61 44 25 100  
 Fax +49 89 61 44 25 200  
 E-Mail: sales@itd-cart.com  
 Web: www.itd-cart.com



### Please note:

Just draw in your selection of accessories at your choice of position between the two dotted lines. Accessories like shelves, drawers etc. for **pro-cart** can be assembled starting **at the first unit (at 32 mm)** from the base.