

vexio-cart Profi

Planning Support

Equipment

Castor diameter:	D _____ mm
Conductivity:	<input type="checkbox"/> yes <input type="checkbox"/> no
Isolating transformers:	<input type="checkbox"/> 660 VA with 4-way EURO socket strip
	<input type="checkbox"/> 1200 VA with 6-way EURO socket strip
Insulating monitor:	<input type="checkbox"/> yes <input type="checkbox"/> no
Socket strip:	ZV. _____
Accessories:	_____ _____ _____

Design

Base colour:	RAL 7035 (Light grey)
Trim colour:	<input type="checkbox"/> .903 (RAL 7035 Light grey)
	<input type="checkbox"/> .908 (RAL 5015 Sky blue)
	<input type="checkbox"/> .909 (RAL 1003 Signal yellow)
	<input type="checkbox"/> _____



35 U (1120 mm)

30 U (960 mm)

25 U (800 mm)

21 U (672 mm)

15 U (480 mm)

10 U (320 mm)

5 U (160 mm)

1 U (32 mm)

Dimensions

Total height TH:

Units (U), 1 U = 32 mm	Castors D 100 mm	Castors D 125 mm
21 U	884 mm	909 mm
30 U	1172 mm	1197 mm
35 U	1332 mm	1357 mm

Total width: 610 mm

Total depth: 620 mm

Project _____

Company _____

Street _____

City _____

Fon _____

Fax _____

E-Mail _____

Processor _____

Please note:

Just draw in your selection of accessories at your choice of position between the two dotted lines. Accessories like shelves, drawers etc. for **pro-cart** can be assembled starting **at the first unit (at 32 mm)** from the base.

ITD GmbH
Lilienthalstr. 19
85579 Neubiberg, Germany
Tel. +49 89 61 44 25 100
Fax +49 89 61 44 25 200
E-Mail: sales@itd-cart.com
Web: www.itd-cart.com

