

Monitor carts

Monitor carts are specially designed to accommodate medical monitors which often need to be moved. When selecting the appropriate working height, integrated mechanical or electrical height adjustment are provided for extra comfort. Ergonomic push handles facilitate manoeuvring and positioning. Four lockable twin castors ensure high comfort while moving and a secure fixture when static.



MONITOR CARTS uni-cart



uni-cart monitor cart with mechanical height adjustment

Equipment cart, 47 U (Total height: 1685 mm), fitted with:

- 1 basic frame with 4 lockable twin castors Ø 100 mm
- 1 mechanical height adjustment
- 1 monitor holder, tiltable with VESA 100 / 200 / 300 adaptation up to 14 kg / 30.8 lbs
- 1 counterweight set
- 2 cable winders, rear
- 1 push handle, rear

Installation height: 47 U (1504 mm)
 Load capacity: 65 kg / 143 lbs
 Finish: RAL 7035 light grey powder-coated
 Support column: anodised
 Decorative trims: push handle see ¹⁾²⁾
 Scope of delivery: assembled

Load capacity of height adjustment	Item number
3-10 kg / 6.6-22 lbs	RS.4914.XXX ¹⁾²⁾
10-16 kg / 22-30.8 lbs	RS.4915.XXX ¹⁾²⁾

uni-cart monitor cart with electric height adjustment

Equipment cart, 35 U (Total height: 1305 mm), fitted with:

- 1 basic frame with 4 lockable twin castors Ø 100 mm
- 1 counterweight set
- 1 electric height adjustment, load capacity up to 20 kg / 44 lbs with power pack and remote control
- 1 monitor holder, universal with VESA 75 / 100 adaptation up to 14 kg / 30.8 lbs
- 1 keyboard drawer with mousepad (left / right)
- 1 EURO socket strip, 3-way with on/off switch
- 2 cable winders, rear
- 1 push handle, rear
- 1 PC / computer holder (small), rear
- 1 connection cable, Schuko, length 5 m
- 1 integrated plug removal barrier

Installation height: 35 U (1120 mm)
 Load capacity: 65 kg / 143 lbs
 Finish: RAL 7035 light grey powder-coated
 Support column: anodised
 Decorative trims: push handle see ¹⁾
 Scope of delivery: assembled

Input voltage	Item number
230 V and 50 / 60 Hz	RS.4145.XXX ¹⁾
110 V and 50 / 60 Hz	RS.4146.XXX ¹⁾

Monitor carts

¹⁾ Decorative trims: please replace .XXX with .901 = RAL 5021 water blue, .902 = RAL 5013 cobalt blue, .903 = RAL 7035 light grey

²⁾ Decorative trims: please replace .XXX with .800 = RAL 7016 anthracite grey

MONITOR CARTS vexio-cart



vexio-cart monitor cart with mechanical height adjustment

Equipment cart, 35 U (Total height: 1320 mm), fitted with:

- 1 basic frame with 4 lockable twin castors Ø 100 mm
- 1 mechanical height adjustment with manual locking
- 1 keyboard drawer with mouse pad (left / right)
- 1 monitor holder, universal with VESA 75/100 adapter, up to 14 kg / 30.8 lbs, and 1 small power pack holder
- 1 push handle, rear
- 2 cable winders, rear
- 1 socket strip, earthed, 3-way with on/off switch
- 2 "Economy" cable ducts, side
- 2 hygiene strips, rear

Installation height: 35 U (1120 mm)
 Load capacity: 65 kg / 143 lbs
 Finish: RAL 7035 light grey powder-coated
 Support column: anodised
 Decorative trims: push handle see ¹⁾
 Scope of delivery: assembled

Load capacity of height adjustment	Item number
3-10 kg / 6.6-22 lbs	NT.5021.XXX ¹⁾
10-16 kg / 22-35.2 lbs	NT.5022.XXX ¹⁾

vexio-cart monitor cart with mechanical height adjustment

Equipment cart, 23 U (Total height: 933 mm), fitted with:

- 1 basic frame with 4 lockable twin castors Ø 100 mm
- 6 hygiene strips
- 1 mechanical height adjustment with release pedal
- 1 monitor holder, tiltable with VESA 100 / 200 / 300 adaptation up to 14 kg / 30.8 lbs with locators
- 1 counterweight set
- 2 cable winders, rear
- 1 push handle, rear

Installation height: 23 U (736 mm)
 Load capacity: 65 kg / 143 lbs
 Finish: RAL 7035 light grey powder-coated
 Support column: anodised
 Decorative trims: push handle / locators see ¹⁾²⁾
 Scope of delivery: assembled

Load capacity of height adjustment	Item number
3-10 kg / 6.6-22 lbs	NT.5014.XXX ¹⁾²⁾
8-14 kg / 17.6-30.8 lbs	NT.5015.XXX ¹⁾²⁾

MONITOR CARTS pro-cart



Monitor carts

pro-cart monitor cart with electric height adjustment

Equipment cart, 35 U (Total height: 1332 mm), fitted with:

- 1 basic frame with 4 locable twin castors Ø 100 mm
- 1 electric height adjustment, load capacity up to 20 kg / 44 lbs with power pack and remote control
- 1 shelf with handle, front
- 1 monitor holder, universal, with VESA 75 / 100 adaptation up to 14 kg / 30.8 lbs
- 1 keyboard drawer with mousepad (left / right)
- 1 push handle, rear
- 2 cable winders, rear
- 1 connection cable, Schuko, length 5 m

Installation height: 35 U (1120 mm)
 Load capacity: 80 kg / 176 lbs
 Finish: RAL 7035 light grey powder-coated
 Decorative trims: handle and push handle see ¹⁾
 Scope of delivery: assembled

Input voltage	Item number
230 V and 50 / 60 Hz	PT.5012.XXX ¹⁾
110 V and 50 / 60 Hz	PT.5017.XXX ¹⁾

pro-cart monitor cart with mechanical height adjustment

Equipment cart, 48 U (Total height: 1798 mm), fitted with:

- 1 basic frame with 4 locable twin castors Ø 125 mm, conductive
- 1 mechanical height adjustment
- 1 monitor holder, tiltable with VESA 100 / 200 / 300 adaptation up to 16 kg / 35.2 lbs
- 1 counterweight set
- 1 push handle, rear
- 2 cable winders, rear

Installation height: 48 U (1536 mm)
 Load capacity: 80 kg / 176 lbs
 Finish: RAL 7035 light grey powder-coated
 Decorative trims: push handle see ¹⁾²⁾
 Scope of delivery: assembled

Load capacity of height adjustment	Item number
3-10 kg / 6.6-22 lbs	PT.5016.XXX ¹⁾²⁾
10-16 kg / 22-35.2 lbs	PT.5015.XXX ¹⁾²⁾

¹⁾ Decorative trims: please replace .XXX with .901 = RAL 5021 water blue, .902 = RAL 5013 cobalt blue, .903 = RAL 7035 light grey

²⁾ Decorative trims: please replace .XXX with .800 = RAL 7016 anthracite grey

MONITOR CARTS FOR PHILIPS IntelliVue



uni-cart rollstand for Philips IntelliVue MP5-70 / MX400-800

Equipment cart, 32 U (Total height: 1265 mm), fitted with:

- 1 basic frame with 4 lockable twin castors Ø 100 mm
- 1 counterweight set
- 1 tilt and swivel unit with Table Top Mount adaptation for Philips IntelliVue MP5-70 / MX400-800 up to 14 kg / 30.8 lbs
- 1 vertical holder for „external rack“
- 1 storage basket
- 1 push handle, rear
- 1 integrated plug removal barrier

Installation height: 32 U (1024 mm)
 Load capacity: 50 kg / 110 lbs
 Finish: RAL 7035 light grey powder-coated
 Decorative trims: push handle and holder: Medium Aqua Accent powder-coated
 Scope of delivery: assembled, Philips Table Top Mount is not included in the scope of delivery

Item number
RS.4903.801

vexio-cart rollstand for Philips IntelliVue MP5-70 / MX400-800

Equipment cart, 25 U (Total height: 1042 mm), fitted with:

- 1 basic frame with 4 lockable twin castors Ø 100 mm
- 1 counterweight set
- 1 tilt and swivel unit with Table Top Mount adaptation for Philips IntelliVue MP5-70 / MX400-800 up to 14 kg / 30.8 lbs
- 1 storage basket
- 1 push handle, rear
- 6 hygiene strips

Installation height: 25 U (800 mm)
 Load capacity: 65 kg / 143 lbs
 Finish: RAL 7035 light grey powder-coated
 Decorative trims: push handle: RAL 7035 light grey powder-coated
 Scope of delivery: assembled, Philips Table Top Mount is not included in the scope of delivery

Item number
NT.5016.903

MONITOR CARTS FOR PHILIPS Efficia



vexio-cart rollstand for Philips Efficia CM 100 / 120 / 150

Equipment cart, 25 U (Total height: 1042 mm),
fitted with:

- 1 basic frame with 4 lockable twin castors Ø 100 mm
- 1 counterweight set
- 1 tilt and swivel unit with Table Top Mount adaptation
for Philips Efficia CM100 / 120 / 150 up to
14 kg / 30.8 lbs
- 1 storage basket
- 1 push handle, rear
- 6 hygiene strips

Installation height: 25 U (800 mm)
Load capacity: 65 kg / 143 lbs
Finish: RAL 7035 light grey powder-coated
Decorative trims: push handle: RAL 7035 light grey
powder-coated
Scope of delivery: assembled, Philips Table Top Mount is
not included in the scope of delivery

Item number

NT.5017.903

