



新一代设备

classic-cart

新产品



技术细节

负载能力

总负载能力	150 kg / 330 lbs
层架的承重能力	50 kg / 110 lbs 分布式负载

尺寸 - 基本框架

宽度 (毫米)	663 / 763 / 1263
深度 (毫米)	660
总高度 (毫米)	840 / 968 / 1128 / 1288 / 1448 / 1608 / 1768

尺寸 - 层架

宽度 (毫米)	500 / 600 / 1100
深度 (毫米)	491

电气化

隔离变压器	1200 VA / 1600 VA / 2000 VA
总开关	yes
插排	yes

接触 >>

Europe: +49 89 614425-0 | sales@itd-cart.com
America: +1 800 947-3901 | salesusa@itd-cart.com
China: +86 21 5031 8818 | saleschina@itd-cart.com
India: +91 86 0015 4064 | salesindia@itd-cart.com
Australia: +61 3 9330 5050 | salesaustralia@itd-cart.com


www.itd-cart.com



定制设计

新颜色

基本色




RAL 9016
traffic white

装饰条颜色



RAL 7016
anthracite grey



RAL 9016
traffic white

新设计

个性化企业设计





DS.0270.799 - 02/2026. 德国印刷。保留所有修改权。

功能优化的高品质设计



新护角 -
设计优雅, 卫生性能更佳



推动手柄 -
手感舒适, 可定制符合
企业品牌标准

平滑的单立柱 -
具有出色的卫生功能
和独特的设计



定制设计 -
识别度高





新的支臂位置 -
活动范围大

额外的背部把手 -
便于从任何位置
对推车进行调整

后门 -
便于存取设备和电缆,
隐蔽式电缆管理

国际插座条 -
畅通无阻



配重 -
隐蔽式安装, 灵活与隔离变压器相结合

